

Case Study: Baxter Healthcare (Thetford)

Pharmaceutical compounding suite for aseptic products

Bassaire successfully constructed the 400 square meter facility on a turnkey contract as well as being the main contractor.

Pharmaceutical grade A, B, C and D areas were validated to MCA standards.

Not only did Bassaire construct the facility but also manufactured the Grade A laminar flow benches using a combination of stainless and powder coated steel.



Other products and services

Bassaire has over fifty years experience in the design, manufacture, installation and maintenance of cleanrooms and clean air products.

